

OrthoLine™ Proximal Femoral Fracture Plating System

Surgical Technique




Arthrex®
Vet Systems

OrthoLine™ Proximal Femoral Fracture Plating System

Introduction


Arthrex has developed many innovations to assist in treating small animal fractures. These new and often unique solutions, combined with state-of-the-art surgeon training, research, and the most well-trained technology specialists in the industry, provide a comprehensive approach for our Arthrex veterinary customers.

The Arthrex OrthoLine proximal femoral fracture (PFF) plating system includes a range of sizes from broad 1.6 mm to broad 3.5 mm. Each plate size is anatomically contoured to mimic the anatomy of patients within a given size range. This feature limits the need for contouring, saving OR time while improving plate placement, bone purchase, and overall screw density. This uniquely designed plate also includes a suture hole to aid in soft tissue closure.




Features and Benefits

- Anatomic shape assists in repair with minimal contouring
- Ideal for subtrochanteric proximal femoral fractures
- Proximal screw trajectories align centrally in metaphyseal and diaphyseal bone
- Anatomic design features include a contour to follow the natural procurvatum of the femur, as well as a dual proximal curvature to follow the trochanter
- Hole spacing and location increase bone purchase options and overall density of screws in the proximal fragment
- Suture hole included on the proximal end of the plate to aid in soft tissue closure
- Locking K-wire guides may be used to estimate screw trajectory prior to drilling
- Scalloped underside to distribute stress and minimize contact
- Multiple K-wire holes and availability of cannulated bending plugs allow for temporary fixation prior to screw insertion

Anatomic Design



Surgical Approach




Use a standard lateral approach to the femur, incise the vastus lateralis at its caudal lateral origin, and reflect cranially.




Reduction

An Intramedullary (IM) pin can be used for initial reduction. The fracture can be further stabilized using the FiberTape® cerclage tension band system. FiberTape cerclage is a strong, simple, reproducible, and effective solution for replacing metal cables and wires traditionally used for fracture management.




Assess Plate

Place the proximal flared portion of the plate against the lateral femur at the origin of the vastus lateralis. Minor contouring through the mid-diaphyseal portion of the plate is often needed.




Contour Plate

If contouring is required, place a threaded bending plug into the locking screw holes where the plate will be contoured. Contour the plate as necessary using bending irons. If large contours are made, take note that screw trajectories may change. If plate length is excessive, the surgeon may cut the plate. Depending on plate size and thickness, this may be performed with either a hand or table-top cutter. If the cutter is not adequate to cut through the plate, a score mark can be made with the cutter, and then bending irons can be used to carefully bend and fatigue the metal until it separates.



Initial Fixation

Once contouring is completed, the bending plugs can be removed or used with K-wires for temporary fixation. BB-Taks and K-wires may also be used in their respective holes.





9

Applying the Plate


Place the first locking screw bicortically in the proximal segment (avoid the reduction IM pin) using the appropriate locking drill guide, drill bit, and depth gauge. Screws may be placed under power and can be brought into contact with the plate. The final turns, however, should be performed manually with the screwdriver.

Note: Surgeon discretion will dictate whether the IM pin is left in place or removed. All screws should be placed to avoid the fracture. The proximal K-wire must be removed to place the third screw from the proximal screw to avoid interference. Final screw placement is at the discretion of the surgeon.






Place the second and third locking screws bicortically in the proximal segment (avoid the reduction IM pin, if present) using the appropriate locking drill guide, drill bit, and depth gauge. It is desirable to capture 6-8 cortices on either side of the fracture.



Place the remaining screws using locking or cortical screws as required.



12

If desired, the suture hole feature of the plate can be used in conjunction with the QuickPass™ SutureLasso™ suture passer to aid in the closure of the vastus lateralis over the plate.

Suture Reference Chart




Plate Size	Plate	VetSuture	Product Description
1.6 mm/2.0 mm		VAR-R316	Polydioxanone 3-0, SH, TP, ½ C
		VAR-R317	Polydioxanone 2-0, SH, TP, ½ C
2.0 mm/2.4 mm		VAR-R334	Polydioxanone 0, CT-2, TP, ½ C
		VAR-R340	Polydioxanone 0, CT-1, TP, ½ C
3.0 mm/3.5 mm		VAR-R468	Polydioxanone 1, CP-1, Rev Ctg, ½ C
		VAR-R468	Polydioxanone 1, CP-1, Rev Ctg, ½ C

Ordering Information

Proximal Femoral Fracture Plates

Product Description	Item Number
1.6 mm Proximal Femoral Broad Plates (Gold)	
Proximal femur plate broad, titanium, 1.6 mm, left (a)	VAR-3116BPF-L
Proximal femur plate broad, titanium, 1.6 mm, right (b)	VAR-3116BPF-R
2.0 mm Proximal Femoral Plates (Blue)	
Proximal femur plate, titanium, 2.0 mm, left (c)	VAR-3120PF-L
Proximal femur plate, titanium, 2.0 mm, right (d)	VAR-3120PF-R
2.4 mm Proximal Femur Plates (Green)	
Proximal femur plate, titanium, 2.4 mm, left (e)	VAR-3124PF-L
Proximal femur plate, titanium, 2.4 mm, right (f)	VAR-3124PF-R
3.0 mm Proximal Femur Plates (Purple)	
Proximal femur plate, titanium, 3.0 mm, left (g)	VAR-3130PF-L
Proximal femur plate, titanium, 3.0 mm, right (h)	VAR-3130PF-R
3.5 mm Proximal Femur Plates (Matte)	
Proximal femur plate, short, SS, 3.5 mm, left (i)	VAR-3035PFS-L
Proximal femur plate, short, SS, 3.5 mm, right (j)	VAR-3035PFS-R
Proximal femur plate, SS, 3.5 mm, left (k)	VAR-3035PF-L
Proximal femur plate, SS, 3.5 mm, right (l)	VAR-3035PF-R
Proximal femur plate broad, SS, 3.5 mm, left (m)	VAR-3035BPF-L
Proximal femur plate broad, SS, 3.5 mm, right (n)	VAR-3035BPF-R



Screws

Product Description	Item Number
1.6 mm Low-Profile Cortical, Variable Angle, Titanium	
Low-profile cortical screw 1.6 mm × 6-20 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20 mm	VAR-8916-06 to -20
Low-profile variable-angle screw 1.6 mm × 6-20 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20 mm	VAR-8916V-06 to -20
2.0 mm Low-Profile Cortical, Locking, Variable Angle, Titanium	
Low-profile cortical screw 2.0 mm × 6-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920-06 to -30
Low-profile locking screw 2.0 mm × 6-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920L-06 to -30
Low-profile variable-angle screw 2.0 mm × 6-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920V-06 to -30
2.4 mm Low-Profile Cortical, Locking, Variable Angle, Titanium	
Low-profile cortical screw 2.4 mm × 8-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924-08 to -30
Low-profile locking screw 2.4 mm × 8-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924L-08 to -30
Low-profile variable-angle screw 2.4 mm × 8-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924V-08 to -30
2.7 mm Low-Profile Cortical, Locking, Stainless Steel	
Low-profile cortical screw 2.7 mm × 10-34 mm Sizes: 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 mm	VAR-8827-10 to -34
Low-profile locking screw 2.7 mm × 10-34 mm Sizes: 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 mm	VAR-8827L-10 to -34
3.0 mm Low-Profile Cortical, Locking, Variable Angle, Titanium	
Low-profile cortical screw 3.0 mm × 8-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930-08 to -40
Low-profile locking screw 3.0 mm × 8-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930L-08 to -40
Low-profile variable-angle screw 3.0 mm × 8-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930V-08 to -40

Screws cont.

Product Description	Item Number
3.5 mm Low-Profile Cortical, Locking, Stainless Steel	
Low-profile cortical screw 3.5 mm × 16-60 mm Sizes: 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60 mm	VAR-8835-16 to -60
Low-profile locking screw 3.5 mm × 16-60 mm Sizes: 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60 mm	VAR-8835L-16 to -60
4.0 mm Low-Profile, Locking, Stainless Steel	
Low-profile locking screw 4.0 mm × 18-60 mm Sizes: 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60 mm	VAR-8840L-18 to -60

Disposables and Limited Reusables


Product Description	Item Number
Drill bit, solid, AO, 1.1 mm (1.6 mm)	VAR-4016D
Drill bit, solid, AO, 1.5 mm (2.0 mm)	VAR-4020D
Drill bit, solid, AO, 1.8 mm (2.4 mm)	VAR-4024D
Drill bit, solid, AO, 2.3 mm (3.0 mm)	VAR-4030D
Drill bit, solid, AO, 2.5 mm (3.5 mm)	VAR-8943-30
Drill bit, solid, AO, 2.8 mm (3.5 mm)	VAR-4035D
Drill bit, solid, AO, 3.5 mm (4.0 mm)	VAR-4040D
Drill bit, solid, short, AO, 1.1 mm (1.6 mm)	VAR-4016SD
Drill bit, solid, short, AO, 1.5 mm (2.0 mm)	VAR-4020SD
Drill bit, solid, short, AO, 1.8 mm (2.4 mm)	VAR-4024SD
Drill bit, solid, short, AO, 2.3 mm (3.0 mm)	VAR-4030SD
Guidewire w/ trocar tip, 0.86 × 80 mm	VAR-8929K
Guidewire w/ trocar tip, 1.1 mm × 150 mm	VAR-8933K
Guidewire w/ trocar tip, 1.3 mm × 150 mm	VAR-8937K

Instruments

Product Description	Item Number
Depth measuring device (1.6 mm/2.0 mm/2.4 mm)	VAR-2024DD
Depth measuring device, (2.7 mm/3.0 mm/3.5 mm/4.0 mm)	VAR-8943-15
T6 driver (1.6 mm/2.0 mm)	VAR-4020-01
T8 driver (2.4 mm)	VAR-4024-01
T10 screwdriver (2.7 mm/3.0 mm)	VAR-8944DH
T15 driver (3.5 mm/4.0 mm)	VAR-8941DH
T6 screwdriver (1.6 mm/2.0 mm)	VAR-4020-02
T8 screwdriver (2.4 mm)	VAR-4024-02
Screw holding forceps (2.7 mm/3.0 mm)	VAR-8943-08
T15 screwdriver (3.5 mm)	VAR-8943-10
Locking plate holder, 2.0 mm	VAR-4020-03
Locking plate holder, 2.4 mm	VAR-4024-03
Locking plate holder, 2.7 mm/3.0 mm	VAR-8950-09

Product Description	Item Number
Locking plate holder, 3.5 mm	VAR-8954-07
Screw holding forceps	VAR-8941F
Drill/depth guide, locking, 1.6 mm	VAR-4016DG
Drill/depth guide, locking, 2.0 mm	VAR-4020DG
Drill/depth guide, locking, 2.4 mm	VAR-4024DG
Drill/depth guide, locking, 2.7 mm	VAR-8950-07
Drill/depth guide, locking, 3.0 mm	VAR-4030DG
Drill/depth guide, locking, 3.5 mm	VAR-4035DG
Drill/depth guide, locking, 4.0 mm	VAR-4040DG
Drill guide, 1.1 mm (1.6 mm)	VAR-4016TDG
Tap/drill guide, 2.0 mm/1.5 mm (2.0 mm)	VAR-4020TDG
Tap/drill guide, 2.4 mm/1.8 mm (2.4 mm)	VAR-4024TDG
2.0 mm/3.0 mm nonlocking drill guide	VAR-8943-31
Tap/drill guide, 3.0 mm/2.3 mm (3.0 mm)	VAR-4030TDG
Drill guide (3.5 mm)	VAR-8943-14
BB-Tak, small, threaded	VAR-8933TBB
BB-Tak, small	VAR-8933BB
BB-Tak, large	VAR-8941BB
BB-Tak, large, threaded	VAR-8941TBB
Drill guide, variable, 1.6 mm	VAR-4016VDG
Drill guide, variable, 2.0 mm	VAR-4020VDG
Drill guide, variable, 2.4 mm	VAR-4024VDG
Drill guide, variable, 3.0 mm	VAR-4030VDG
Bone tap, 2.0 mm	VAR-4020T
Bone tap, 2.4 mm	VAR-4024T
Bone tap, 3.0 mm	VAR-4030T
K-wire drill guide, 0.86 mm (1.6 mm/2.0 mm)	VAR-4020KDG
K-wire drill guide, 1.14 mm (2.4 mm)	VAR-4024KDG
K-wire drill guide, 1.14 mm (2.7 mm/3.0 mm)	VAR-4030KDG
K-wire drill guide, 1.3 mm (3.5 mm)	VAR-4035KDG
Bending plug, cannulated, 1.6 mm/2.0 mm	VAR-4020-04
Bending plug, cannulated, 2.4 mm	VAR-4024-04
Bending plug, cannulated, 3.0 mm	VAR-4030-04
Bending plug, cannulated, 3.5 mm	VAR-4035-04
Bending iron, small (1.6 mm/2.0 mm)	VAR-4000-07
Bending iron, medium (2.4 mm/3.0 mm)	VAR-4000-08
Bending iron, large (3.5 mm/3.5 mm broad)	VAR-4000-09
Freer elevator	VAR-4000-10
Hohmann retractor, double ended, 6 mm/10 mm	VAR-4000-11
Ikuta clamp	VAR-4000-12
Lobster clamp, mini	VAR-4000-13
Lobster clamp, mini, radiolucent	VAR-4000-14
Periosteal elevator, 6 mm curved blade	VAR-4000-15
Pliers, needlenose	VAR-4000-16
Pointed reduction forceps	VAR-4000-17
Reduction forceps, guidewire	VAR-4000-18
Sharp hook	VAR-4000-19
Termite forceps	VAR-4000-20
Toothed reduction forceps, Kocher	VAR-4000-21

Cases and Caddies

Image	Product Description	Item Number
 <p>The image shows a rectangular, silver-colored metal case with a perforated top and bottom. The top surface features the 'Arthrex' logo and 'VET ORTHOLINE™ SYSTEM' text. The bottom surface has 'eDFU' and 'VAR-4000GC' printed on it.</p>	<p>OrthoLine™ system case</p>	<p>VAR-4000GC</p>
 <p>The image shows a white plastic insert with a grid of holes. Various metal surgical instruments are mounted on the insert. Labels include 'VAR-4000GC-01', 'eDFU', 'SCREWDRIVER', 'DRIVER', 'WIRE BENDER', 'WIRE CUTTER', 'BALL DRIVER', 'LOCKING BALL DRIVER', 'BALL END', and 'BALL END'. The Arthrex logo is also visible.</p>	<p>Generic case insert</p>	<p>VAR-4000GC-01</p>
 <p>The image shows a black plastic caddy with a grid of holes. A grey metal tray with the Arthrex logo and 'eDFU' is shown next to it. The tray has a grid of holes and a handle.</p>	<p>1.6 mm Screw caddy</p>	<p>VAR-3016SC-01</p>
 <p>The image shows a white plastic caddy with a grid of holes. A grey metal tray with the Arthrex logo and 'eDFU' is shown next to it. The tray has a grid of holes and a handle.</p>	<p>2.0 mm Screw caddy</p>	<p>VAR-3020SC-01</p>
 <p>The image shows a black plastic caddy with a grid of holes. A grey metal tray with the Arthrex logo and 'eDFU' is shown next to it. The tray has a grid of holes and a handle.</p>	<p>2.4 mm Screw caddy</p>	<p>VAR-3024SC-01</p>

Cases and Caddies

Image	Product Description	Item Number
 <p>The image shows a black plastic screw caddy with a grid of holes. A clear plastic lid is placed next to it, featuring the Arthrex logo and 'eDFU' branding.</p>	<p>3.0 mm Screw caddy</p>	<p>VAR-3030SC-01</p>
 <p>The image shows a black plastic screw caddy with a grid of holes. A clear plastic lid is placed next to it, featuring the Arthrex logo and 'eDFU' branding.</p>	<p>3.5 mm/4.0 mm Screw caddy</p>	<p>VAR-4035SC-02</p>
 <p>The image shows a black plastic bending plug caddy with a grid of holes. A clear plastic lid is placed next to it, featuring the Arthrex logo.</p>	<p>Bending plug caddy</p>	<p>VAR-4000BPC</p>



This is not veterinary advice and Arthrex recommends that veterinarians be trained in the use of any particular product before using it in surgery. A veterinarian must always rely on his or her own professional clinical judgment when deciding whether to use a particular product. A veterinarian must always refer to the package insert, product label and / or instructions for use before using any Arthrex product. Products may not be available in all markets because product availability is subject to the regulatory or veterinary practices in individual markets. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. Please contact your Arthrex representative if you have questions about availability of products in your area.

arthrexvetsystems.com

© 2024-06 Arthrex, Inc. All rights reserved. vLT1-000297-en-US_A



Arthrex manufacturer,
authorized representative,
and importer information
(Arthrex eIFUs)



US patent information