

OrthoLine™ Distal Humeral Fracture System

Surgical Technique



Arthrex®
Vet Systems

OrthoLine™ Distal Humeral Fracture System


Introduction

The first anatomic medial distal humeral fracture plate to address medial, supracondylar, T, and Y fractures with a single plate is available in sizes ranging from 2.0 mm to 3.5 mm. OrthoLine plates were thoughtfully engineered by incorporating surgeon feedback and designed for a range of patients.

Features and Benefits

- Plate incorporates the transcondylar screw
- Optimized trajectory for the transcondylar screw
- 4 screws in the distal aspect
- Anatomic plate design with left and right options
- Strong single-bridging plate design
- Epicondylar relief
- Designed for medial, supracondylar, T, and Y fractures¹


Anatomic Design





Use a standard bilateral approach to the distal humerus.

Note: The medial approach may incorporate the modified Wendelburg surgical approach and is available in video format on the Arthrex Vet Systems website. Reduce the fracture and hold with a temporary reduction pin placed through the point of the epicondyle.



Place the plate on the medial surface of the bone with the point of the epicondyle and the temporary reduction pin (if used) placed over and through the large central aperture. It is important that the plate is generally placed as far caudal and proximal as the epicondyle will allow, which should facilitate alignment

of the transcondylar screw. Minimal contouring may be required to fit the anatomy, but **this may alter desired screw trajectory** and thus should be kept to a minimum. During contouring, place cannulated bending plug inserts into the locking screw holes where applicable.




The plate may be temporarily affixed to the bone using a K-wire, cannulated bending plug, or threaded BB-Tak (the BB-Tak should only be used on the proximal aspect). It is important to note the distal K-wire hole has the same trajectory as the transcondylar screw and thus may be used as a reference guide. To achieve appropriate trajectory, hand place the K-wire to rest in the hole prior to driving with the pin driver.




4

Place the K-wire through the locking K-wire guide and confirm the appropriate trajectory by palpating the pin or with a fluoroscopic image. If the trajectory is acceptable, use the appropriate drill bit to create the pilot hole for locking screw insertion. Prior to screw insertion, place a Vusellum or point-to-point forceps to aid in fracture compression.


If the K-wire placement is not acceptable, redirect the placement and use the drill bit to drill a pilot hole for a nonlocking cortical screw or, alternatively, on titanium plates (3.0 mm and lower), use a variable angle-locking screw.




4a




4b



4c






4d



Place the third proximal screw. First, lock in the locking drill guide, drill using the appropriate drill bit, measure, and place the screw. Screws may be placed under power. The final turns, however, should be performed manually with the screwdriver.

Note: The variable-angle guide can also be used for variable-angle locking titanium screws sizes 3.0 mm and below.

Drill Steps: Set 2





Start at the most distal screw on the distal end of the plate and work proximally. For each screw, lock the locking drill guide in the screw hole, drill using the appropriate drill, measure, and place the screw. Screws may be placed under power. The final turns, however, should be performed manually with the manual screwdriver. If the fracture site is below a screw hole, do not place a screw in this location, a bending plug can be used to fill the hole. K-wires can be removed at this point.


Note: The variable-angle guide can also be used for variable-angle locking titanium screws sizes 3.0 mm and below.

Note: The screws in this step are downsized. For example, the 3.5 mm plate uses 2.7 mm screws, the 3.0 mm plate uses 2.4 mm screws, the 2.4 mm plate uses 2.0 mm screws, and the 2.0 mm plate uses either 2.0 mm or 1.6 mm screws.



Apply the proximal screws. Start at the most distal screw on the proximal side of the fracture line and work proximal. First, lock in the locking drill guide, drill using the appropriate drill bit, measure, and place the screw. Screws may be placed under power. The final turns, however, should be performed manually with the manual screwdriver. Remove the BB-Tak and place the next locking screw.

Note: The variable-angle guide can also be used for variable-angle locking titanium screws sizes 3.0 mm and below.



Caudal view




Medial view

Ordering Information

Distal Humeral Plates

Product Description	Item Number
Distal Humeral Plate, Ti, 2.0 mm, left (a)	VAR-3120DH-L
Distal Humeral Plate, Ti, 2.0 mm, right (b)	VAR-3120DH-R
Distal Humeral Plate, Ti, 2.4 mm, left (c)	VAR-3124DH-L
Distal Humeral Plate, Ti, 2.4 mm, right (d)	VAR-3124DH-R
Distal Humeral Plate, Ti, 3.0 mm, left (e)	VAR-3130DH-L
Distal Humeral Plate, Ti, 3.0 mm, right (f)	VAR-3130DH-R
Distal Humeral Plate, SS, 3.5 mm, left (g)	VAR-3035DH-L
Distal Humeral Plate, SS, 3.5 mm, right (h)	VAR-3035DH-R



Screws

Product Description	Item Number
1.6 mm Low-Profile Cortical, Variable Angle	
Low-Profile Cortical Screw, 1.6 mm × 6 mm-20 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20 mm	VAR-8916-06 to -20
Low-Profile Variable-Angle Screw, 1.6 mm × 6 mm-20 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20 mm	VAR-8916V-06 to -20
2.0 mm Low-Profile Cortical, Variable Angle, Locking	
Low-Profile Cortical Screw, 2.0 mm × 6 mm-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920-06 to -30
Low-Profile Variable-Angle Screw, 2.0 mm × 6 mm-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920V-06 to -30
Low-Profile Locking Screw, 2.0 mm × 6 mm-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920L-06 to -30
2.4 mm Low-Profile Cortical, Variable Angle, Locking	
Low-Profile Cortical Screw, 2.4 mm × 8 mm-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924-08 to -30
Low-Profile Variable-Angle Screw, 2.4 mm × 8 mm-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924V-08 to -30
Low-Profile Locking Screw, 2.4 mm × 8 mm-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924L-08 to -30
2.7 mm Low-Profile Cortical, Locking	
Low-Profile Cortical Screw, 2.7 mm × 10 mm-34 mm Sizes: 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 mm	VAR-8827-10 to -34
Low-Profile Locking Screw, 2.7 mm × 10 mm-34 mm Sizes: 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 mm	VAR-8827L-10 to -34
3.0 mm Low-Profile Cortical, Variable Angle, Locking	
Low-Profile Cortical Screw, 3.0 mm × 8 mm-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930-08 to -40
Low-Profile Variable-Angle Screw, 3.0 mm × 8 mm-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930V-08 to -40
Low-Profile Locking Screw, 3.0 mm × 8 mm-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38 40 mm	VAR-8930L-08 to -40

Screws cont.

Product Description	Item Number
3.5 mm Low-Profile Cortical, Locking	
Low-Profile Cortical Screw, 3.5 mm × 8 mm-65 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65 mm	VAR- 8835-08 to - 65
Low-Profile Locking Screw, 3.5 mm × 10 mm-60 mm Sizes: 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60 mm	VAR- 8835L-10 to - 60
4.0 mm Low-Profile Cortical, Locking	
Low-Profile Locking Screw, 4.0 mm × 18 mm-65 mm Sizes: 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60, 65 mm	VAR- 8840L-18 to - 65

Equipment Part List

Product Description	Item Number
Disposable/Limited Reusable	
Drill Bit, solid, AO, 1.1 mm	VAR- 4016D
Drill Bit, solid, AO, 1.5 mm	VAR- 4020D
Drill Bit, solid, AO, 1.8 mm	VAR- 4024D
Drill Bit, solid, AO, 2.3 mm	VAR- 4030D
Drill Bit, solid, short, AO, 2.5 mm	VAR- 8943-30
Drill Bit, solid, 2.7 mm	VAR- 8944-22
Drill Bit, solid, AO, 2.8 mm	VAR- 4035D
Drill Bit, solid, AO, 3.5 mm	VAR- 4040D
Drill Bit, solid, short, AO, 1.1 mm	VAR- 4016SD
Drill Bit, solid, short, AO, 1.5 mm	VAR- 4020SD
Drill Bit, solid, short, AO, 1.8 mm	VAR- 4024SD
Drill Bit, solid, short, AO, 2.3 mm	VAR- 4030SD
Guidewire w/ Trocar, 0.86 mm × 80 mm	VAR- 8929K
Guidewire w/ Trocar, 1.1 mm × 150 mm	VAR- 8933K
Guidewire w/ Trocar, 1.3 mm × 150 mm	VAR- 8937K
Instruments	
Depth Measuring Device, 1.6 mm/2.0 mm/2.4 mm	VAR- 2024DD
Depth Measuring Device, 2.7 mm/3.0 mm/3.5 mm/4.0 mm	VAR- 8943-15
T6 Driver, 1.6 mm/2.0 mm	VAR- 4020-01
T8 Driver, 2.4 mm	VAR- 4024-01
T10 Driver, 3.0 mm	VAR- 8944DH
T15 Driver, 3.5 mm	VAR- 8941DH
T6 Screwdriver, 1.6 mm/2.0 mm	VAR- 4020-02
T8 Screwdriver, 2.4 mm	VAR- 4024-02
T10 Screwdriver, 2.7 mm/3.0 mm	VAR- 8943-08
T15 Screwdriver, 3.5 mm/4.0 mm	VAR- 8943-10
Locking Plate Holder, 2.0 mm	VAR- 4020-03
Locking Plate Holder, 2.4 mm	VAR- 4024-03
Locking Plate Holder, 2.7 mm/3.0 mm	VAR- 8950-09
Locking Plate Holder, 3.5 mm	VAR- 8954-07
Screw Holding Forceps	VAR- 8941F

Product Description	Item Number
Drill/Depth Guide, locking, 1.6 mm	VAR- 4016DG
Drill/Depth Guide, locking, 2.0 mm	VAR- 4020DG
Drill/Depth Guide, locking, 2.4 mm	VAR- 4024DG
Drill/Depth Guide, locking, 3.0 mm	VAR- 4030DG
Drill/Depth Guide, locking, 3.5 mm	VAR- 4035DG
Drill/Depth Guide, locking, 4.0 mm	VAR- 4040DG
Drill Guide, 1.1 mm	VAR- 4016TDG
Tap/Drill Guide, 1.5 mm	VAR- 4020TDG
Tap/Drill Guide, 1.8 mm	VAR- 4024TDG
Tap/Drill Guide, 2.0 mm	VAR- 8943-31
Tap/Drill Guide, 2.3 mm	VAR- 4030TDG
Tap/Drill Guide, 2.4 mm	VAR- 8943-14
BB-Tak, small, threaded	VAR- 8933TBB
BB-Tak, small	VAR- 8933BB
BB-Tak, large, threaded	VAR- 8941TBB
BB-Tak, large	VAR- 8941BB
Drill Guide, variable, 1.6 mm	VAR- 4016VDG
Drill Guide, variable, 2.0 mm	VAR- 4020VDG
Drill Guide, variable, 2.4 mm	VAR- 4024VDG
Drill Guide, variable, 3.0 mm	VAR- 4030VDG
Bone Tap, 2.0 mm	VAR- 4020T
Bone Tap, 2.4 mm	VAR- 4024T
Bone Tap, 2.7 mm	VAR- 4027T
Bone Tap, 3.0 mm	VAR- 4030T
K-Wire Drill Guide, 0.86 mm (1.6 mm/2.0 mm)	VAR- 4020KDG
K-Wire Drill Guide, 1.14 mm (2.4 mm)	VAR- 4024KDG
K-Wire Drill Guide, 1.14 mm (2.7 mm/3.0 mm)	VAR- 4030KDG
K-Wire Drill Guide, 1.3 mm (3.5 mm)	VAR- 4035KDG
Bending Plug, cannulated, 1.6/2.0 mm	VAR- 4020-04
Bending Plug, cannulated, 2.4 mm	VAR- 4024-04
Bending Plug, cannulated, 2.7 mm	VAR- 4027-04
Bending Plug, cannulated, 3.0 mm	VAR- 4030-04
Bending Plug, cannulated, 3.5 mm	VAR- 4035-04
Bending Iron, small, 1.6 mm/2.0 mm	VAR- 4000-07
Bending Iron, medium, 2.4 mm/3.0 mm	VAR- 4000-08
Bending Iron, large 3.5 mm/broad 3.5 mm	VAR- 4000-09
Freer Elevator	VAR- 4000-10
Hohmann Retractor, double ended, 6 mm/10 mm	VAR- 4000-11
Ikuta Clamp	VAR- 4000-12
Lobster Clamp, mini	VAR- 4000-13
Lobster Clamp, mini, radiolucent	VAR- 4000-14
Periosteal Elevator, 6 mm curved blade	VAR- 4000-15
Pliers, needle nose	VAR- 4000-16
Pointed Reduction Forceps	VAR- 4000-17
Reduction Forceps, guidewire	VAR- 4000-18
Sharp Hook	VAR- 4000-19
Termite Forceps	VAR- 4000-20
Toothed Reduction Forceps, Kocher	VAR- 4000-21

Cases and Caddies

Image	Product Description	Item Number
	OrthoLine™ Case	VAR-4000GC
	Generic Case Insert	VAR-4000GC-01
	1.6 mm Screw Caddy	VAR-3016SC-01
	2.0 mm Screw Caddy	VAR-3020SC-01
	2.4 mm Screw Caddy	VAR-3024SC-01
	2.7 mm Screw Caddy	VAR-4027SC-01

Cases and Caddies cont.

Image	Product Description	Item Number
	3.0 mm Screw Caddy	VAR-3030SC-01
	3.5 mm/4.0 mm Screw Caddy	VAR-4035SC-02
	Bending Plug Caddy	VAR-4000BPC

Products advertised in this surgical technique guide may not be available in all countries. For information on availability, please contact Arthrex Customer Service or your local Arthrex representative.

Reference

1. Arthrex, Inc. Data on file (APT 05162). Naples, FL; 2021.



This is not veterinary advice and Arthrex recommends that veterinarians be trained in the use of any particular product before using it in surgery. A veterinarian must always rely on their own professional clinical judgment when deciding whether to use a particular product. A veterinarian must always refer to the package insert, product label, and/or directions for use before using any Arthrex product. Products may not be available in all markets because product availability is subject to the regulatory or veterinary practices in individual markets. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. Please contact your Arthrex representative if you have questions about availability of products in your area.

View U.S. patent information at www.arthrex.com/corporate/virtual-patent-marking

arthrexvetsystems.com