

OrthoLine™ Ilial Fracture System

Surgical Technique




Arthrex®
Vet Systems

OrthoLine™ Ilial Fracture System


Introduction

The OrthoLine ilial fracture plate system includes a range of sizes from 1.6 mm to 3.5 mm, with an additional long 3.5 plate for larger dogs. Each plate size is anatomically contoured to mimic the anatomy of patients within a given size range. The plate includes cranial divergent screw trajectories to improve purchase in the cranial ilium. The caudal screw trajectories are more dorsally oriented to help diverge away from the coxofemoral joint. Additionally, the ilial plate includes a suture hole that fits Arthrex VetSuture for apposition of the gluteal musculature at closure.

Features and Benefits

- Divergent cranial screw trajectories for added pull-out strength
- Caudal screw trajectories to diverge from coxofemoral joint
- Designed for cranial, caudal, and oblique fracture patterns
- Anatomic plate design with left and right options
- High screw density where needed
- Suture hole to aid in soft tissue closure
- Scalloped underside to distribute stress and minimize contact

Anatomic Design




Fracture Patterns




Transverse fracture pattern



Cranial fracture pattern




Caudal fracture pattern



Long oblique fracture pattern


Surgical Technique



1

To determine contour requirements, place the implant on the bone surface of the ilium. The plate may be shifted cranial or caudal based on surgeon preference and specific fracture type. In plate application, the most-ventral screws can be placed just cranial to the cranial


pillar of the acetabulum. A caudal twist may be required at the most-caudal aspect of the plate. Care must be taken with large plate contours as screw trajectories will be altered.



2


Using the appropriate screwdriver, place the cannulated threaded bending plugs into the locking screw holes where the plate will be contoured. Contour the plate as necessary using bending irons or another form of plate benders.

To bend the tip of the plate, use the fork end of the bending iron. Once the contouring is complete, the bending plugs can be removed or used in conjunction with the appropriate size K-wire for temporary fixation.




3

The plate may be temporarily affixed to the bone using multiple methods. For temporary fixation, 3 methods may be used: K-wire holes, K-wires in conjunction with cannulated bending plugs, and/or the application of a threaded BB-Tak in the center of a universal hole.




First, affix the plate to the cranial and caudal fracture segments using a cortical or locking screw. When using the locking guide, it is important to note that the guide must be inserted parallel to the screw hole. This is not always perpendicular to the plate surface in Arthrex anatomic plates, as screw trajectories are

unique. Inappropriate drill guide alignment may result in locking thread damage and should be avoided. Using a second sterile plate as a reference guide in drilling, measuring, and placing screws.




Note: If desired, the variable-angle guide can be used for titanium implants sizes 3.0 mm and below. If care is not taken, a screw placed using a variable-angle method can interfere with another screw.



Place locking screws caudal to the fracture as required by drilling using the locking drill guide, measuring, and placing the screw. When using power, the screw is brought into contact with the plate and should be manually locked into place. Ensure the screws avoid the fracture line.


Note: The variable-angle guide can also be used for variable-angle locking titanium screws sizes 3.0 mm and below.




Place locking screws cranial of the fracture as required by drilling using the locking drill guide, measuring, and placing the screw. The screw is brought into contact with the plate and should be manually locked into place. Screw trajectories may be varied to incorporate the sacral body for maximum fixation strength of the cranial segment. Titanium variable-angle locking titanium screws may be used in sizes 3.0 mm and below, or a

cortical screw can be placed in any of the universal holes in the Arthrex Iliac fracture plate. Ensure the screws avoid the fracture line.

Note: The variable-angle guide can also be used for variable-angle locking titanium screws sizes 3.0 mm and below.




Secure VetSuture to the gluteal muscle in a mattress pattern, cut the needle from the suture, feed the suture end through the suture hole, and tie over the plate. A QuickPass™ SutureLasso™ suture passer may be used to aid in placement of the VetSuture. Finish with routine closure. The QuickPass SutureLasso suture passer is placed under the suture hole in the ilial plate. The nitinol loop is advanced by rolling the thumbwheel from top to bottom. Ensure the SutureLasso suture passer




has passed through the suture hole, then place the VetSuture into the lasso and move the thumbwheel forward to tension the suture passer. Pull the suture through the suture hole and tie it in a routine fashion.


Note: If the plate is placed against the bone, the suture may need to be placed prior to bringing the plate to the bone with screw fixation.




Finish with routine closure.



Lateral view.




Medial screw trajectory view, with outline of the sacral body.



11

Caudal view.












12

Depending on the fracture pattern, placement of the plate can be cranial or caudal.

Suture Reference Chart



Plate Size	Plate	VetSuture	Product Description
1.6 mm/2.0 mm		VAR-R316	Polydioxanone 3-0, SH, TP, 1/2 C
		VAR-R317	Polydioxanone 2-0, SH, TP, 1/2 C
		VAR-J8665	Polypropylene 3-0, FS-2, Rev Ctg, 3/8 C
2.0 mm/2.4 mm		VAR-R334	Polydioxanone 0, CT-2, TP, 1/2 C
		VAR-R340	Polydioxanone 0, CT-1, TP, 1/2 C
		VAR-R467	Polydioxanone 0, CP-1, Rev Ctg, 1/2 C
3.0 mm/3.5 mm		VAR-R468	Polydioxanone 1, CP-1, Rev Ctg, 1/2 C
			
			

Ordering Information

Iliac Fracture Plates

Product Description	Item Number
1.6 mm Iliac Fracture Plates (Gold)	
Iliac fracture plate, titanium, 1.6 mm, left (a)	VAR-3116I-L
Iliac fracture plate, titanium, 1.6 mm, right (b)	VAR-3116I-R
2.0 mm Iliac Fracture Plates (Blue)	
Iliac fracture plate, titanium, 2.0 mm, left (c)	VAR-3120I-L
Iliac fracture plate, titanium, 2.0 mm, right (d)	VAR-3120I-R
2.4 mm Iliac Fracture Plates (Green)	
Iliac fracture plate, titanium, 2.4 mm, left (e)	VAR-3124I-L
Iliac fracture plate, titanium, 2.4 mm, right (f)	VAR-3124I-R
3.0 mm Iliac Fracture Plates (Purple)	
Iliac fracture plate, titanium, 3.0 mm, left (g)	VAR-3130I-L
Iliac fracture plate, titanium, 3.0 mm, right (h)	VAR-3130I-R
3.5 mm Iliac Fracture Plates (Matte)	
Iliac fracture plate, SS, 3.5 mm, left (i)	VAR-3035I-L
Iliac fracture plate, SS, 3.5 mm, right (j)	VAR-3035I-R
Iliac fracture plate, long, SS, 3.5 mm, left (k)	VAR-3035IL-L
Iliac fracture plate, long, 3.5 mm, right (l)	VAR-3035IL-R



(a)



(b)



(c)



(d)



(e)



(f)



(g)



(h)



(i)



(j)



(k)



(l)

Screws

Product Description	Item Number
1.6 mm Low-Profile Cortical, Variable Angle, Titanium	
Low-profile cortical screw 1.6 mm × 6-20 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20 mm	VAR-8916-06 to -20
Low-profile variable-angle screw 1.6 mm × 6-20 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20 mm	VAR-8916V-06 to -20
2.0 mm Low-Profile Cortical, Locking, Variable Angle, Titanium	
Low-profile cortical screw 2.0 mm × 6-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920-06 to -30
Low-profile locking screw 2.0 mm × 6-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920L-06 to -30
Low-profile variable-angle screw 2.0 mm × 6-30 mm Sizes: 6, 7, 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8920V-06 to -30
2.4 mm Low-Profile Cortical, Locking, Variable Angle, Titanium	
Low-profile cortical screw 2.4 mm × 8-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924-08 to -30
Low-profile locking screw 2.4 mm × 8-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924L-08 to -30
Low-profile variable-angle screw 2.4 mm × 8-30 mm Sizes: 8, 9, 10, 11, 12, 13, 14, 16, 18, 20, 22, 24, 26, 28, 30 mm	VAR-8924V-08 to -30
2.7 mm Low-Profile Cortical, Locking, Stainless Steel	
Low-profile cortical screw 2.7 mm × 10-34 mm Sizes: 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 mm	VAR-8827-10 to -34
Low-profile locking screw 2.7 mm × 10-34 mm Sizes: 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34 mm	VAR-8827L-10 to -34
3.0 mm Low-Profile Cortical, Locking, Variable Angle, Titanium	
Low-profile cortical screw 3.0 mm × 8-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930-08 to -40
Low-profile locking screw 3.0 mm × 8-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930L-08 to -40
Low-profile variable-angle screw 3.0 mm × 8-40 mm Sizes: 8, 10, 12, 14, 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40 mm	VAR-8930V-08 to -40

Screws cont.

Product Description	Item Number
3.5 mm Low-Profile Cortical, Locking, Stainless Steel	
Low-profile cortical screw 3.5 mm × 16-60 mm Sizes: 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60 mm	VAR-8835-16 to -60
Low-profile locking screw 3.5 mm × 16-60 mm Sizes: 16, 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60 mm	VAR-8835L-16 to -60
4.0 mm Low-Profile, Locking, Stainless Steel	
Low-profile locking screw 4.0 mm × 18-60 mm Sizes: 18, 20, 22, 24, 26, 28, 30, 32, 34, 36, 38, 40, 42, 44, 46, 48, 50, 52, 54, 56, 58, 60 mm	VAR-8840L-18 to -60

Disposables and Limited Reusables

Product Description	Item Number
Drill bit, solid, AO, 1.1 mm (1.6 mm)	VAR-4016D
Drill bit, solid, AO, 1.5 mm (2.0 mm)	VAR-4020D
Drill bit, solid, AO, 1.8 mm (2.4 mm)	VAR-4024D
Drill bit, solid, AO 2.0 (2.7 mm)	VAR-8944-22
Drill bit, solid, AO, 2.3 mm (3.0 mm)	VAR-4030D
Drill bit, solid, AO, 2.5 mm (3.5 mm)	VAR-8943-30
Drill bit, solid, AO, 2.8 mm (3.5 mm)	VAR-4035D
Drill bit, solid, AO, 3.5 mm (4.0 mm)	VAR-4040D
Drill bit, solid, short, AO, 1.1 mm (1.6 mm)	VAR-4016SD
Drill bit, solid, short, AO, 1.5 mm (2.0 mm)	VAR-4020SD
Drill bit, solid, short, AO, 1.8 mm (2.4 mm)	VAR-4024SD
Drill bit, solid, short, AO, 2.3 mm (3.0 mm)	VAR-4030SD
Guidewire w/ trocar tip, 0.86 mm × 80 mm	VAR-8929K
Guidewire w/ trocar tip, 1.1 mm × 150 mm	VAR-8933K
Guidewire w/ trocar tip, 1.3 mm × 150 mm	VAR-8937K

Instruments

Product Description	Item Number
Depth measuring device (1.6 mm/2.0 mm/2.4 mm)	VAR-2024DD
Depth measuring device (2.7 mm/3.0 mm/3.5 mm/4.0 mm)	VAR-8943-15
T6 driver (1.6 mm/2.0 mm)	VAR-4020-01
T8 driver (2.4 mm)	VAR-4024-01
T10 screwdriver (2.7 mm/3.0 mm)	VAR-8944DH
T15 driver (3.5 mm/4.0 mm)	VAR-8941DH
T6 screwdriver (1.6 mm/2.0 mm)	VAR-4020-02
T8 screwdriver (2.4 mm)	VAR-4024-02
Screw holding forceps (2.7 mm/3.0 mm)	VAR-8943-08
T15 screwdriver (3.5 mm)	VAR-8943-10
Locking plate holder, 2.0 mm	VAR-4020-03
Locking plate holder, 2.4 mm	VAR-4024-03
Locking plate holder, 2.7 mm/3.0 mm	VAR-8950-09
Locking plate holder, 3.5 mm	VAR-8954-07
Screw holding forceps	VAR-8941F

Product Description	Item Number
Drill/depth guide, locking, 1.6 mm	VAR-4016DG
Drill/depth guide, locking, 2.0 mm	VAR-4020DG
Drill/depth guide, locking, 2.4 mm	VAR-4024DG
Drill/depth guide, locking, 2.7 mm	VAR-8950-07
Drill/depth guide, locking, 3.0 mm	VAR-4030DG
Drill/depth guide, locking, 3.5 mm	VAR-4035DG
Drill/depth guide, locking, 4.0 mm	VAR-4040DG
Drill guide, 1.1 mm (1.6 mm)	VAR-4016TDG
Tap/drill guide, 2.0 mm/1.5 mm (2.0 mm)	VAR-4020TDG
Tap/drill guide, 2.4 mm/1.8 mm (2.4 mm)	VAR-4024TDG
2.0 mm/3.0 mm nonlocking drill guide	VAR-8943-31
Tap/drill guide, 3.0 mm/2.3 mm (3.0 mm)	VAR-4030TDG
Drill guide (3.5 mm)	VAR-8943-14
BB-Tak, small, threaded	VAR-8933TBB
BB-Tak, small	VAR-8933BB
BB-Tak, large	VAR-8941BB
BB-Tak, large, threaded	VAR-8941TBB
Drill guide, variable, 1.6 mm	VAR-4016VDG
Drill guide, variable, 2.0 mm	VAR-4020VDG
Drill guide, variable, 2.4 mm	VAR-4024VDG
Drill guide, variable, 3.0 mm	VAR-4030VDG
Bone tap, 2.0 mm	VAR-4020T
Bone tap, 2.4 mm	VAR-4024T
Bone tap, 2.7 mm	VAR-4027T
Bone tap, 3.0 mm	VAR-4030T
K-wire drill guide, 0.86 mm (1.6 mm/2.0 mm)	VAR-4020KDG
K-wire drill guide, 1.14 mm (2.4 mm)	VAR-4024KDG
K-wire drill guide, 1.14 mm (2.7 mm/3.0 mm)	VAR-4030KDG
K-wire drill guide, 1.3 mm (3.5 mm)	VAR-4035KDG
Bending plug, cannulated, 1.6 mm/2.0 mm	VAR-4020-04
Bending plug, cannulated, 2.4 mm	VAR-4024-04
Bending plug, cannulated, 2.7 mm	VAR-4027-04
Bending plug, cannulated, 3.0 mm	VAR-4030-04
Bending plug, cannulated, 3.5 mm	VAR-4035-04
Bending iron, small (1.6 mm/2.0 mm)	VAR-4000-07
Bending iron, medium (2.4 mm/3.0 mm)	VAR-4000-08
Bending iron, large (3.5 mm/3.5 mm broad)	VAR-4000-09
Freer elevator	VAR-4000-10
Hohmann retractor, double ended, 6 mm/10 mm	VAR-4000-11
Ikuta clamp	VAR-4000-12
Lobster clamp, mini	VAR-4000-13
Lobster clamp, mini, radiolucent	VAR-4000-14
Periosteal elevator, 6 mm curved blade	VAR-4000-15
Pliers, needlenose	VAR-4000-16
Pointed reduction forceps	VAR-4000-17
Reduction forceps, guidewire	VAR-4000-18
Sharp hook	VAR-4000-19
Termite forceps	VAR-4000-20
Toothed reduction forceps, Kocher	VAR-4000-21

Cases and Caddies

Image	Product Description	Item Number
 <p>The image shows a rectangular, silver-colored metal case with a perforated top and bottom. The top surface features the 'Arthrex' logo and 'VET ORTHOLINE SYSTEM' text. The bottom surface also has the 'Arthrex' logo and 'VET ORTHOLINE SYSTEM' text, along with 'VAR-4000GC' and 'eDFU' branding.</p>	OrthoLine™ case	VAR-4000GC
 <p>The image shows a white, rectangular insert tray with a grid of holes. It is labeled 'VAR-4000GC-01' and 'eDFU'. Various orthopedic tools are mounted in the holes, including screws, ball sockets, and other components. Labels like 'SCREWDRIVER', 'DRIVER', 'HANDLE', 'BALL SOCKET', and 'LOCKING BALL' are visible.</p>	Generic case insert	VAR-4000GC-01
 <p>The image shows a black, rectangular screw caddy with a handle. It is designed to hold 1.6 mm screws. The caddy is shown next to its corresponding grey insert tray, which has a grid of holes for the screws. The 'Arthrex' and 'eDFU' logos are visible on the insert tray.</p>	1.6 mm Screw caddy	VAR-3016SC-01
 <p>The image shows a white, rectangular screw caddy with a handle. It is designed to hold 2.0 mm screws. The caddy is shown next to its corresponding grey insert tray, which has a grid of holes for the screws. The 'Arthrex' and 'eDFU' logos are visible on the insert tray.</p>	2.0 mm Screw caddy	VAR-3020SC-01
 <p>The image shows a black, rectangular screw caddy with a handle. It is designed to hold 2.4 mm screws. The caddy is shown next to its corresponding grey insert tray, which has a grid of holes for the screws. The 'Arthrex' and 'eDFU' logos are visible on the insert tray.</p>	2.4 mm Screw caddy	VAR-3024SC-01

Cases and Caddies

Image	Product Description	Item Number
 <p>A black plastic screw caddy with a grid of holes for screws. The lid is open, showing a white insert with a grid of holes. The Arthrex logo and 'RDFU' are visible on the lid.</p>	3.0 mm Screw caddy	VAR-3030SC-01
 <p>A white plastic screw caddy with a grid of holes for screws. The lid is open, showing a white insert with a grid of holes. The Arthrex logo and 'RDFU' are visible on the lid.</p>	2.7 mm Screw caddy	VAR-4027SC-01
 <p>A black plastic screw caddy with a grid of holes for screws. The lid is open, showing a white insert with a grid of holes. The Arthrex logo and 'RDFU' are visible on the lid.</p>	3.0 mm Screw caddy	VAR-3030SC-01
 <p>A black plastic screw caddy with a grid of holes for screws. The lid is open, showing a white insert with a grid of holes. The Arthrex logo and 'RDFU' are visible on the lid.</p>	3.5 mm/4.0 mm Screw caddy	VAR-4035SC-02
 <p>A black plastic bending plug caddy with a grid of holes for screws. The lid is open, showing a white insert with a grid of holes. The Arthrex logo and 'RDFU' are visible on the lid.</p>	Bending plug caddy	VAR-4000BPC



This is not veterinary advice and Arthrex recommends that veterinarians be trained in the use of any particular product before using it in surgery. A veterinarian must always rely on his or her own professional clinical judgment when deciding whether to use a particular product. A veterinarian must always refer to the package insert, product label and / or instructions for use before using any Arthrex product. Products may not be available in all markets because product availability is subject to the regulatory or veterinary practices in individual markets. Postoperative management is patient-specific and dependent on the treating professional's assessment. Individual results will vary and not all patients will experience the same postoperative activity level or outcomes. Please contact your Arthrex representative if you have questions about availability of products in your area.

arthrexvetsystems.com

© 2024-06 Arthrex, Inc. All rights reserved. vLT1-000301-en-US_A



Arthrex manufacturer,
authorized representative,
and importer information
(Arthrex eIFUs)



US patent information