


3.5 mm TPLO System


VAR-4035S



Pic.	Item Number	Qty.	Description
1	VAR-4035D		Drill Bit, 2.8 mm, long
2	VAR-4000-CAL		TPLO Caliper
3	VAR-8943-14		Nonlocking Drill Guide, 2.5 mm/3.5 mm
4	VAR-8943-30		Drill Bit, 2.5 mm, short
5	VAR-8943-15		Depth Device, low profile, 2.7 mm/3.5 mm
6	VAR-8943-10		T15 Screwdriver
7	VAR-8941DH		T15 Driver
8	VAR-4035DG		Locking Drill Guide, 3.5 mm
9	VAR-4035SC-01		TPLO Screw Caddy, 3.5 mm
10	VAR-8941F		Screw Holding Forceps
11	VAR-8835-16 to -50		Low-Profile Cortical Screw, 3.5 mm × 16 mm to 50 mm
12	VAR-8835L-16 to -50		Low-Profile Locking Screw, 3.5 mm × 16 mm to 50 mm

3.5 mm TPLO System

VAR-4035S



Pic.	Item Number	Qty.	Description
14	VAR-8967K		K-Wire, .092 in (2.3 mm)
15	VAR-13303-3.0		Guide Pin, .118 in (3.0 mm)
16	VAR-8933K		K-Wire, .045 in (1.14 mm)
17	VAR-8941K		K-Wire, .062 in (1.57 mm)
18	VAR-4000TPSM		TPLO Tissue Protector, small
19	VAR-8954-07		Locking Plate Holder, 3.5 mm
20	VAR-4000TPS		TPLO Tissue Protector, standard
21	VAR-4035SM-L and -R		TPLO Locking Plate, 3.5 mm, small, left and right
22	VAR-4035S-L and -R		TPLO Locking Plate, 3.5 mm, standard, left and right
23	VAR-4035B-L and -R		TPLO Locking Plate, 3.5 mm, broad, left and right

3.5 mm TPL0 System

VAR-4035S

Date of Surgery:

Surgeon Name:

Account:

Representative Contact Information:

Set Total:

Miscellaneous Notes:

To obtain a copy of a product's directions for use (DFU), which includes sterilization instructions, please visit our website at www.arthrex.com/edfu or contact your local representative.

Page 3 of 3

www.arthrexvetsystems.com

© 2019 Arthrex, Inc. All rights reserved. vDOC1-000143-en-US_A



Arthrex[®]
CANINE
InternalBrace™ SYSTEM
Vet Systems