











Distal Radial Fracture Plate Size Chart





OrthoLine™ System Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Radial Plate, Ti, 1.6 mm		VAR-3116DR-L	39.0	2.6	1.4	7.3	5.0
		VAR-3116DR-R					
Distal Radial Plate, Ti, Broad, 1.6 mm		VAR-3116BDR-L	42.3	3.1	1.8	7.9	5.4
		VAR-3116BDR-R					
Distal Radial Plate, Ti, 2.0 mm		VAR-3120DR-L	46.9	3.7	1.9	8.8	5.7
		VAR-3120DR-R					
Distal Radial Plate, Ti, 2.4 mm		VAR-3124DR-L	57.2	4.7	2.5	10.0	7.1
		VAR-3124DR-R					



Distal Radial Fracture Plate Size Chart

OrthoLine™ System Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Radial Plate, Ti, 3.0 mm		VAR-3130DR-L	64.1	5.4	3.2	11.6	8.8
		VAR-3130DR-R					
Distal Radial Plate, SS, 3.5 mm		VAR-3035DR-L	80.1	6.3	4.1	13.6	11.4
		VAR-3035DR-R					