











# OrthoLine™ System Fracture Plate





## Size Chart for Iliac Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Iliac Fracture Plate, Ti, 1.6 mm		VAR-3116I-L	35.8	1.8	2.0	7.9	5.0
		VAR-3116I-R					
Iliac Fracture Plate, Ti, 2.0 mm		VAR-3120I-L	38.0	2.5	2.7	9.1	6.1
		VAR-3120I-R					
Iliac Fracture Plate, Ti, 2.4 mm		VAR-3124I-L	45.3	3.3	3.5	11.6	7.6
		VAR-3124I-R					
Iliac Fracture Plate, Ti, 3.0 mm		VAR-3130I-L	60.8	3.9	5.2	13.4	9.4
		VAR-3130I-R					



# OrthoLine™ System Fracture Plate

## Size Chart for Iliac Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Iliac Fracture Plate, SS, 3.5 mm		VAR-3035I-L	71.0	4.3	6.2	14.7	10.9
		VAR-3035I-R					
Iliac Fracture Plate, Long, SS, 3.5 mm		VAR-3035IL-L	94.3	4.3	6.4	14.7	10.9
		VAR-3035IL-R					