









OrthoLine™ TPLO Plates

Size Chart

	Product Image	Part Number	Length (mm)	Thickness (mm)	Max Height (mm)
TPLO Plate, Ti, 1.6 mm		VAR-4116S-L	24.6	1.9	6.1
		VAR-4116S-R			
TPLO Plate, Ti, 2.0 mm		VAR-4120S-L	29.0	2.3	6.1
		VAR-4120S-R			
TPLO Plate, Ti, 2.4 mm		VAR-4124S-L	35.4	2.7	7.8
		VAR-4124S-R			
TPLO Plate, SS, 2.7 mm		VAR-4027S-L	43.4	2.9	8.3
		VAR-4027S-R			
TPLO Plate, Ti, 3.0 mm		VAR-4030S-L	48.8	3.3	8.3
		VAR-4030S-R			

OrthoLine™ TPLO Plates

Size Chart

	Product Image	Part Number	Length (mm)	Thickness (mm)	Max Height (mm)
TPLO Plate, SS, Small, 3.5 mm		VAR-4035SM-L	53.5	3.7	9.7
		VAR-4035SM-R			
TPLO Plate, SS, Standard, 3.5 mm		VAR-4035S-L	62.4	3.7	10.7
		VAR-4035S-R			
TPLO Plate, SS, Broad, 3.5 mm		VAR-4035B-L	76.6	3.7	11.6
		VAR-4035B-R			
TPLO Plate, SS, 4.5 mm		VAR-4045S-L	89.2	3.7	12.2
		VAR-4045S-R			

