

OrthoLine Fracture Management System

Weight Reference Chart



| | 1-4 kg | 4-8 kg | 8-12 kg | 13-19 kg | 20-30 kg | 31-45 kg | 45-55+ kg |
|----------------|-----------------|----------------|----------------|----------------|----------|-----------------|-----------|
| Radius | 1.6 mm, 1.6b mm | 2.0 mm | 2.0 mm, 2.4 mm | 3.0 mm | 3.5 mm | 3.5 mm | 3.5 mm |
| Humerus | 1.6 mm, 1.6b mm | 2.0 mm | 2.4 mm, 3.0 mm | 3.0 mm, 3.5 mm | 3.5 mm | 3.5b mm | 3.5b mm |
| Femur | 1.6 mm, 1.6b mm | 2.0 mm, 2.4 mm | 2.4 mm, 3.0 mm | 3.0 mm, 3.5 mm | 3.5 mm | 3.5 mm, 3.5b mm | 3.5b mm |
| Tibia | 1.6 mm, 1.6b mm | 2.0 mm | 2.4 mm | 3.0 mm, 3.5 mm | 3.5 mm | 3.5 mm, 3.5b mm | 3.5b mm |
| T-Plate | 1.6 mm, 1.6b mm | 2.0 mm | 2.0 mm, 2.4 mm | 3.0 mm | 3.5 mm | 3.5 mm | 3.5 mm |

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.