



3.0 mm TPLO System

VAR-4030S



Pic.	Item Number	Qty.	Description
1	VAR-8967K		K-Wire, 2.35 mm (0.092 in)
2	VAR-8937K		K-Wire, 1.3 mm (0.051 in)
3	VAR-8933K		K-Wire, 1.1 mm (0.045 in)
4	VAR-13203-05		Osteotome, 5 mm
5	VAR-8826M		Mini mallet, 7 in
6	VAR-4000TPSM		Tissue protector, small
7	VAR-13210		Hohmann retractor
8	VAR-8950-09		Plate holder, 2.7 mm/3.0 mm
9	VAR-4130S-R		TPLO plate, 3.0 mm standard, right
10	VAR-4130S-L		TPLO plate, 3.0 mm standard, left





Pic.	Item Number	Qty.	Description
16	VAR-4030DG		Locking drill guide, 3.0 mm
17	VAR-4030TDG		Nonlocking drill guide, 3.0 mm
18	VAR-4000-CAL		TPLO caliper
19	VAR-4030VDG		Drill guide, variable, 3.0 mm
20	VAR-4030D		Drill bit, solid, AO, 2.3 mm
21	VAR-4030SD		Drill bit short, solid, AO, 2.3 mm
22	VAR-8943-15		Depth device
23	VAR-8943-08		T10 screwdriver
24	VAR-8944DH		T10 driver
25	VAR-8930-08 to -40		Low-profile cortical screw, titanium, 3.0 mm × 8 mm to 40 mm
26	VAR-8930V-08 to -40		Low-profile variable angle screw, titanium, 3.0 mm × 8 mm to 40 mm
27	VAR-8930L-08 to -40		Low-profile locking screw, titanium, 3.0 mm × 8 mm to 40 mm
28	VAR-8941F		Screw holding forceps
29	VAR-3030SC-01		3.0 mm screw caddy



3.0 mm TPLO System

VAR-4030S

Date of Surgery:

Surgeon Name:

Account:

Representative Contact Information:

Set Total:

Miscellaneous Notes:

To obtain a copy of a product's directions for use (DFU), which includes sterilization instructions, please visit our website at www.arthrex.com/edfu or contact your local representative.



Arthrex®
CANINE
InternalBrace™ SYSTEM

Arthrex®
Vet Systems