




Drill Torque Speed Weight

Reference Chart

Product Image	Drill	Handpiece Weight	Drill Speed	Drill Torque	Reamer Speed	Reamer Torque	Cutting Speed
	300	550 g (1.2 lb)	1300 rpm	24 Ncm (2 in lb)	500 rpm	1.15 Nm (10 in lb)	15,000 rpm
	400	630 g (1.3 lb)	1300 rpm	3.2 Nm (28 in lb)	330 rpm	9 Nm (79 in lb)	22,000 rpm
	600	950 g (2.09 lb)	950 rpm	4.5 Nm (40 in lb)	200 rpm	17 Nm (150 in lb)	13,000 rpm