









Distal Humeral Plate Size Chart

OrthoLine™ System Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Humeral Plate Ti, 2.0 mm		VAR-3120DH-L	48.2	4.2	2.4	12.7	6.0
		VAR-3120DH-R					
Distal Humeral Plate Ti, 2.4 mm		VAR-3124DH-L	55.8	5.6	3.2	15.2	7.6
		VAR-3124DH-R					
Distal Humeral Plate Ti, 3.0 mm		VAR-3130DH-L	75.5	8.3	4.0	20.3	9.1
		VAR-3130DH-R					
Distal Humeral Plate SS, 3.5 mm		VAR-3035DH-L	85.2	9.2	4.4	22.3	10.6
		VAR-3035DH-R					

