











Distal Femur Fracture Plate Size Chart



OrthoLine™ System Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Plate, Ti, Broad, 1.6 mm		VAR-3116BDF-L	68.3	1.9	1.8	6.1	5.0
		VAR-3116BDF-R					
Distal Femur Plate, Ti, 2.0 mm		VAR-3120DF-L	81.6	3.1	2.3	7.4	6.0
		VAR-3120DF-R					
Distal Femur Plate, Ti, 2.4 mm		VAR-3124DF-L	104.4	3.5	2.8	9.6	7.3
		VAR-3124DF-R					
Distal Femur Plate, Ti, 3.0 mm		VAR-3130DF-R	131.8	3.6	3.6	10.6	9.3
		VAR-3130DF-L					



Distal Femur Fracture Plate Size Chart

OrthoLine™ System Fracture Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Distal Femur Plate, SS, 3.5 mm		VAR-3035DF-L	161.3	4.2	4.3	13.05	11.9
		VAR-3035DF-R					