


Small and Mighty

Mini PushLock® CCL Anchor for Stifle Stabilization



- The smallest CCL anchor on the market
- 2.5 mm diameter and 8 mm length designed for patients under 10 lb (4.5 kg)
- Radiolucent PEEK construct
- Reliable strength in knotless and knotted techniques¹
- Available in convenience pack

Lateral Suture Stabilization



See our surgical
technique video
and guide at:

ArthrexVetSystems.com

Reference

1. Roca Hernandez R. Ex vivo mechanical properties of a 2.5 mm bone anchor for treatment of cranial cruciate ligament rupture in toy breed dogs [unpublished data]. 2020.



Arthrex®
Vet Systems