

# TPLO Plate System

## Weight Reference Chart



Internal/Brace™ Technique Anchors	TPLO Plate Sizes									Internal/Brace Technique Sutures
2.5 mm PushLock® anchor	1.6 mm									#2 FiberWire® suture SutureTape suture, 1.3 mm
2.4 mm/2.5 mm PushLock anchor		2.0 mm								SutureTape suture, 1.3 mm
2.5 mm/2.9 mm PushLock anchor			2.4 mm							SutureTape suture, 1.3 mm
3.5 mm SwiveLock® anchor				2.7 mm						LabralTape™ suture, 1.5 mm FiberTape® suture, 2.0 mm
3.5 mm/4.75 mm SwiveLock anchor					3.0 mm					LabralTape suture, 1.5 mm FiberTape suture, 2.0 mm
3.5 mm/4.75 mm SwiveLock anchor						3.5 mm, small				LabralTape suture, 1.5 mm FiberTape suture, 2.0 mm
4.75 mm SwiveLock anchor							3.5 mm, standard			FiberTape suture, 2.0 mm
4.75 mm/5.5 mm SwiveLock anchor								3.5 mm, broad		FiberTape suture, 2.0 mm
5.5 mm SwiveLock anchor									4.5 mm	FiberTape suture, 2.0 mm

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultant's recommendation. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.