


# Canine *Internal/Brace*™ Instrument Case


VAR-1679C



Pic.	Item Number	Description
1	VAR-1679C	Outer Case With Pin Mat
2	VAR-1679C-01	TightRope® Tray
3	VAR-1679C-02	Anchor Tray
4	VAR-1679C-03	Meniscal Tray

# Canine *Internal/Brace*™ Instrument Case

TightRope® Tray, VAR-1679C-01



Pic.	Item Number	Description	Mini TightRope® CCL	Standard TightRope CCL	Mini TightRope Hip	Standard Hip
1	VAR-1529	Tensioner	✓	✓	✓	✓
2	VAR-11797	Small Joint FiberWire® Scissors	✓	✓	✓	✓
3	VAR-8911DC	2.7 mm Drill Bit, cannulated (optional TPLO tibial tunnel)	✓		✓	
4	VAR-8920DC	3.5 mm Drill Bit, cannulated (for double 5.5 mm)		✓		✓
5	VAR-8920P	K-Wire, 1.24 mm (0.049 in)	✓	✓	✓	✓
6	VAR-8935G	Drill Guide, 3.5 mm/4.0 mm		✓		✓
7	VAR-2810	Aiming Guide	✓	✓	✓	✓


## TightRope Implant Weight Sizes

Mini CCL	15-40 lbs (6-18 kg)
Standard CCL	35-250 lbs (15-113 kg)
Mini Hip	20-50 lbs (9-22 kg)
Standard Hip	50+ lbs (22+ kg)



# Canine *Internal/Brace*™ Instrument Case


Anchor Tray, VAR-1679C-02



Pic.	Item Number	Description	2.5 mm PushLock® Anchor	3.5 mm SwiveLock® Anchor	4.75 mm SwiveLock Anchor	5.5 mm SwiveLock Anchor
1	VAR-2324PTB	Tap, 4.75 mm			✓	
2	VAR-1927CTB	Tap, 5.5 mm				✓
3	VAR-1678-03	Tap, 3.5 mm		✓		
4	VAR-2923D-01	PushLock Drill Bit, 2.9 mm (coming soon)	N/A	N/A	N/A	N/A
5	VAR-2923D-02	PushLock Drill Bit Hard Bone, 2.9 mm (coming soon)	N/A	N/A	N/A	N/A
6	VAR-1927D	Spade Tip Drill			✓	✓
7	VAR-1678-05	3.5 mm SwiveLock Drill Bit		✓		
8	VAR-8737-09	2.5 mm Drill Bit, cannulated		✓	✓	✓
9	VAR-8933-20C	2.0 mm Drill Bit, cannulated	✓			
10	VAR-2922D-01	2.4 PushLock Drill Bit (coming soon)	N/A	N/A	N/A	N/A
11	VAR-8920DC	3.5 mm Drill Bit, cannulated (for double 5.5 mm)				✓
12	VAR-8911DC	2.7 mm Drill Bit, cannulated (optional TPLO tibial tunnel)			✓	✓
13	VAR-11797	Small Joint FiberWire® Scissors	✓	✓	✓	✓
14	VAR-2924DG	2.9/2.4 mm PushLock Anchor Drill Guide				
15	VAR-1678-01	SwiveLock Drill Guide, 3.5 mm		✓		
16	VAR-8933K	K-Wire, 1.14 mm (0.045 in)	✓	✓	✓	✓
17	VAR-1255-08	Suture Passing Wire	✓	✓	✓	✓
18	VAR-8941K	K-Wire, 1.57 mm (0.062 in)				

# Canine *Internal/Brace*™ Instrument Case

Meniscal Tray, VAR-1679C-03



Pic.	Item Number	Description
1	VAR-4000	Meniscal Blade Cutter
2	VAR-5003	Push Blade
3	VAR-5005	Hook Blade
4	VAR-4010-27	Stifle Distractor, large
5	VAR-4010-20	Stifle Distractor, medium
6	VAR-4010-13	Stifle Distractor, small
7	VAR-5020	Cannulated Handle
8	VAR-30000	Small Joint Probe
9	VAR-4001	Meniscal Blade Cutter, instruments

## Stifle Distractor Weight Sizes

Small	20-40 lbs (9-18 kg)
Medium	35-75 lbs (16-34 kg)
Large	75-120 lbs (34-54 kg)