











# OrthoLine™ System Fracture Plate







## Size Chart for Proximal Femur Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Proximal Femur Plate, Ti, Broad, 1.6 mm		VAR-3116BPF-L	64.2	7.0	1.9	7.0	5.0
		VAR-3116BPF-R					
Proximal Femur Plate, Ti, 2.0 mm		VAR-3120PF-L	78.6	7.9	2.2	8.1	6.1
		VAR-3120PF-R					
Proximal Femur Plate, Ti, 2.4 mm		VAR-3124PF-L	88.8	9.7	2.9	9.5	7.6
		VAR-3124PF-R					
Proximal Femur Plate, Ti, 3.0 mm		VAR-3130PF-L	154.3	11.5	3.6	12.6	9.4
		VAR-3130PF-R					



# OrthoLine™ System Fracture Plate

## Size Chart for Proximal Femur Plate

	Product Image	Part Number	Length (mm)	Max Thickness Head (mm)	Max Thickness Shaft (mm)	Max Width Head (mm)	Max Width Shaft (mm)
Proximal Femur Plate, SS, Short, 3.5 mm		VAR-3035PFS-L	155.1	13.0	4.3	15.2	11.9
		VAR-3035PFS-R					
Proximal Femur Plate, SS, 3.5 mm		VAR-3035PF-L	188.1	13.0	4.3	15.2	11.9
		VAR-3035PF-R					
Proximal Femur Plate, SS, Broad, 3.5 mm		VAR-3035BPF-L	201.7	14.3	5.5	17.7	14.0
		VAR-3035BPF-R					