

















MPL Plate System

Weight Reference Chart

								
	1-4 kg	3-7 kg	8-12 kg	13-17 kg	15-20 kg	18-25 kg	26-39 kg	40-55 kg
Plate Size	1.6 mm	2.0 mm	2.4 mm	2.7 mm	3.0 mm	3.5 mm, small	3.5 mm, standard	3.5 mm, broad
Plate								

All weight ranges are suggestions based on the ultimate load of the product and/or on our veterinary consultants' recommendations. It is up to the veterinarian to determine the proper product, application, technique, and weight range for using the selected product.